

BT-A95 Electro Hydraulic Operating Table (Side -controlled)



It is for surgeries of the general surgery department, heart renal, orthopedics, neurosurgery, gynecology and urology. The tabletop lifting is provided by a hydraulic pump, and all the adjustments of the table are on the both sides.

The tabletop and outer cover can use carbon steel with plastic spraying or high quality stainless steel, based on your needs.

Dimension: 2010*500*800-1050mm

Main Function:

☐ Lifting range: 800-1050mm

☐ Lean forward: 35°

☐ Lean backward: 22°

☐ Lean left: 22°

☐ Lean right: 22°

☐ Back board folded up: 75°

☐ Back board folded down: 22°

☐ Head board folding up: 45°

☐ Head board folding down: 90°

Leg board folding down: 90°

Main Feature:

1. Single layer table surface or double optional
2. The rise and drop of operating table was handled by pedaling the oil pump, other various movements were handled by the mechanical transmission head.
3. The leg board is detachable, abducted and folded, easy to adjust, and very convenient for various operations;
4. The head and leg plates are detachable for easy operation and provide more flexible operation space.
5. Memory Cotton Pad:Antibacterial and easy to clean,waterproof grade Ip8.
6. Standardized side rail :Universal block locking for all kinds of infusion stand, anesthesia Stand.
7. The whole is made of 304 stainless steel ,smooth shape,the main table and removable sub table for special position surgery.and the surface is treated with special matt, which is more modern and easier to clean.
8. Weight capacity:250 kg
9. Warranty period:1 year

Standard Accessory	
Leg support	1 pair
Arm supports	1 pair
Shoulder and waist holders	2 pairs
Anesthesia shelf	1 pcs



Dedication Beyond Measure

+254 720 941 445 / +254 789 671 699

+254 722 525 079 / +254 733 859 265

villasurgicals@yahoo.com | villasurgicals2030@gmail.com

P.O. Box 13163 - 00100, Nairobi, Kenya

Viken Thirty Industrial Park, Along Eastern Bypass, Kenya