



RO-CHAIN



STAR OF OT LIGHT SUPPLIER
LED EXAMINATION SHADOWLESS LIGHT

A light used during examinations or treatment of the patient. It is usually used during examination or treatment during minor events
Designed to be easily moved from one location to another.

Arm: 105 +/- cm articulated, spring loaded arm, arm with on/off switch and incorporated electronically transformer.

Medical Light to be proposed as basic equipment in health structures, it can be used for medical and gynaecological examination and minor operation.

One of the advantages of a LED examination lamp is its compactness. It can be placed in front of the operator without obstructing sight.



SHANGHAI RO-CHAIN MEDICAL CO.,LTD

Address: Building 1, Kangqiao Rd 787, Pudong New District, Shanghai, China

Tel: 0086-21-58055329 Fax: 0086-21-58055330


Mobile(WhatsApp):+0086-13248293519

E-mail: SevenY@ro-chainmed.com

Official website: www.ro-chainmed.com



RC-YD01A(LED)



Main Technical Data

RC-YD01A(LED)	LED Examination lamp
Illuminance	30,000Lux 0.5m/5,000Lux 1m
Colour temperature	4800±200K
Adjust	adjust light
Total irradiance	380W/ m²
Size of light field	100mm
Rated power of bulb	3.3/3W
Power supply voltage	220V/50Hz
bulb	LED bulb

Main Technical Data

RC-YD01A	Halogen Examination lamp
Illuminance	20W:20,000Lux 0.5m/5,000Lux 1m
Colour temperature	3800±200K
Adjust	adjust light
Total irradiance	380W/ m²
Size of light field	100mm
Rated power of bulb	12V/20W
Power supply voltage	220V/50Hz
bulb	halogen bulb