

Global Markets



T6 | Turbine-driven
ICU & Transport Ventilator

An intelligent ventilation solution for all patients and situations

Compliance with the standards of emergency and transport ventilator

and intensive care ventilator

All-In-One Ventilator for intensive care units and emergency transport

- A minimum tidal volume of 2ml
- Work for all patient groups like adults, pediatrics, and neonates
- Invasive and non-invasive ventilation



Pre-hospital and
Intra-hospital Transport



ICU Treatment



Pediatric and
Neonatal Ventilation



Small yet Powerful

- Lightweight at less than 6.2kg, easy to carry with just one hand
- 10.4-inch color touchscreen, clear and intuitive
- Allow for a visual display of patient

data on a secondary monitor, separate from the primary screen
· Rugged and durable, ideal choice for extreme

environments like shocks, vibrations, rain, and dust

Turbine-driven with Optimal Performance

Advanced turbine technology

- A powerful, quiet turbine that offers constant ventilation

Quick, intuitive use

- Clear and easy-to-use interface, simplifying your work and saving time to therapy

Safe and reliable measurement

- Proximal flow sensor solution allows to measure flow and pressure close to the patient's airway with high accuracy



All the conventional and advanced ventilation modes give you a variety of therapy options in one device

Invasive and Non-invasive Ventilation


- Volume controlled ventilation modes: IPPV, V-A/C, V-SIMV
- Pressure controlled ventilation modes: PCV, P-A/C, P-SIMV,
- Spontaneous mode: CPAP/PSV
- Advanced modes: PRVC, PRVC-SIMV, APRV, BiPPV
- High-flow oxygen therapy mode: HFNC

Configurations	•Standard •Optional
Aminimumtidalvolumeof2ml	◦
BasicVentilationModes: V-A/C,IPPV,V-SIMV,PCV,P-A/C, P-SIMV,CPAP/PSV,CPR	•
AdvancedVentilationModes: PRVC,PRVC-SIMV,APRV,BIPPV	◦
HFNC	◦
Trends,waveforms,loops,CPR,suctioning maneuver,medicationnebulization,dynamiclung	•
P-Vtool,lungrecruitment	◦
Transportandstationarysolution (Trolley,humidifier,supportarm)	◦
EtCO2module	◦

Note: Configurations in the table above are subject to change without prior notice.



Continuous real-time monitoring of patient during the entire ventilation period





Intuitive, easy-to-use and versatile

Ergonomic and Intuitive User Interface

-Allow easy and precise control of waveforms, loops, and parameters settings

Advanced Functions

-HFNC, CPR, lung recruitment, dynamic lung, P-V tool, and P0.1

Innovative and Safe Treatment

-The system is completed with a medical device AC power supply and internal high-capacity rechargeable battery. It also allows operation and rapid charging with the power source of a ground vehicle

-High-pressure and low-pressure oxygen sources

Excellent After Sales Support

-Detachable breathing valve to prevent cross-infection after autoclave sterilization

-Medical flow sensor allows to be disinfected many times

-Bedside hook supports for easy transfer

Integrated Trolley Support System

- Sleek, stable and highly mobile, this trolley is an all-in-one, comprehensive solution
- Easy to assemble and disassemble components on the trolley
- Hold a full range of compatible accessories
- Provide great convenience and ease-of-use when ventilating and moving patients

Easy to install

The trolley is designed to easily hold or remove ventilator and accessories, including oxygen cylinders, a humidifier, a circuit support arm and breathing circuits.



CO2 module



Humidifier



Two emergency oxygen cylinders



Breathing circuits

